


# MAINTENANCE OF YOUR WALKER.

## Brake Adjusting

1. To ensure proper braking, from time to time adjust the brake pads located on the top of the rear wheels. To adjust the brake, remove the rear facing wheel cover, or the Slow Down Brake if installed, by removing the 2 Phillips screws on front of the rear fork.
2. Using a 2.5 mm hex key, loosen the set screw and move the brake pad down for a higher braking force. Test the brakes making sure that when applying the parking brakes, the rear wheels stay in the lock position.
3. Reinstall the rear fork cover or Slow Down Brake if installed.



## Handle Height Adjustment

To adjust the handles to a suitable height, stand behind the walker with your shoulders level. The handles should be at the same height or slightly higher than the user's wrists. A therapist might suggest a different handle height to accommodate a specific need for the user.

To adjust the handles, remove the knob and pull the hex bolt out. Adjust the handle to the prescribed or desired handle height and insert the hex bolt ensuring that the hexagonal part of the bolt fits properly on the hexagonal cut out on the

**Caution: Do not over tighten.**

## Installing Basket

To install the folding soft basket, fold the Xpresso part way by pulling on the center seat handle and insert the basket black brackets on the rails on both inside frames, push on the middle of the seat to open the walker. **Caution: Ensure you hearing the "Click" sound for proper positioned the basket.**

