

P E R S O N A L T R A N S P O R T C H A I R



Seat Colors
Black Vinyl

Body Color
Red

Roll

**Disassembles
Quickly &
Easily!**

GOLDEN
LiteRider[®]ptc

The LiteRider PTC is easier than ever to transport! This new design creates a lightweight base that disassembles into two pieces with the heaviest piece weighing only 35 lbs! Model GP162 offers a compact turning radius with a robust 300 lb. weight capacity. For the consumer looking for high quality at a great value, nothing beats the all-new LiteRider PTC!



GOLDEN
technologies

YOUR LIFE IN MOTION

Model Number	GP162	
Weight Capacity	300 lbs.	
Drive Wheels	Rear	
Maximum Speed	3.5 mph	
Operating Range ¹	9.5 miles	
Ground Clearance	2.5"	
Turning Radius	26"	
Type Batteries	2-18 AH	
Colors	Red	
Freewheel Mode	Yes	
Electronic Speed Control	Yes	
Electro-Mechanical Brakes	Yes	
Charger	Off Board	
Controller	PG VR-2 50 Amp Remote	
Dimensions:		
Length	33"	
Width	22.5"	
Height (ground to top of back of seat)	33" - 34"	
Height (ground to top of back of headrest)	N/A	
Ground to Top of Seat	20.5" - 21.5"	
Number of Seat Height Positions	2	
Size of Increments	1.0"	
Weight of Unit: (assembled)	116 lbs.	
Front Half	34 lbs.	
Rear Half	35 lbs.	
Seat w/Arms	22 lbs.	
Batteries	25 lbs/pack	
Tires:	Flat Free	
Drive Wheels	9"	
Casters	6"	
Anti-Tip Wheels	3"	
Standard Seat:	Stadium Style	
Slide Seat	N/A	
Back Height - no headrest	16"	
Back Height with headrest	N/A	
Width x Depth	17" x 16"	
Color	Black Vinyl	
Optional Seat:	Stadium Style	Stadium Style
Back Height - no headrest	16"	16"
Width x Depth	18" x 16"	20" x 18"
Color	Black Vinyl	Black Vinyl
Width x Depth	17" x 16"	
Footrest:	Yes	
Height adjustable	Yes	
Angle adjustable	Yes	
Warranty:		
Drive Train	13 Month Warranty	
Electronics ²	13 Month Warranty	
Frame	Five Years	

¹ Battery range will vary due to rider weight, drive surface and terrain.
² Electronics warranty excludes batteries.

Disassembles easily for transport!



1. Remove seat



2. Pull out battery pack



3. Grasp the frame and locking handle; lift the handle to unlock, and lift the frame up and off the drive train.



Perfect for transport!

Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Touch-up Paint
- Deluxe Pack N' Go
- Seat Belt (Optional)

Some accessories require a mounting bracket or mounting clips.