



## **AVAILABLE FABRIC**





## **AVAILABLE TRIM COLORS**



Radiant Red

Sunburst

Orange



Power Pink

Galactic

Green



Ocean Blue



Midnight Black

## PREMIUM TRIM COLORS



Mossy Oak® Breakup® Country™



Carbon Fiber

## **SPECIFICATIONS**

| Maximum User Weight        | 300 lbs  |
|----------------------------|--|
| Maximum Speed              | 6 mph<br>5 mph <i>without anti-tippers</i>                         |
| Range                      | 13 mi*   |
| Base Width                 | 24"  |
| Base Length                | 40" with anti-tippers 36" without anti-tippers                     |
| Minimum Turning Radius     | 26.5"  |
| Ground Clearance           | 3"   |
| Recommended Battery Type   | Group 34 & Group 24  |
| Weight Including Batteries | 386 lbs with Group 34 batteries<br>405 lbs with Group 24 batteries |
| Drive Electronics          | R-Net 120A   |

| Seat-to-Floor Height (STFH) | 17.5 , 18.5 , 19.5  |
|-----------------------------|---|
| Optional Seat Elevator      | 12"   |
| Tilt Options                | 0° - 50° <i>Posterior</i><br>5°, 10°, 20° & 30° <i>Anterior</i> |
| Recline Options             | 85° - 120° (manual)<br>85° - 180° (power)                       |
| Legrest Angle               | 85° - 170°  |
| Backrest Height             | 20", 23" - 28" (by 1" incr.)                                    |
| Seat Depth                  | 14" - 22" (by 1" incr.)   |
| Seat Width                  | 17", 19", 21", 23"  |
| Armrest Pad Lengths         | 10", 13", 16", 18"  |
|                             |   |

 $<sup>{}^{\</sup>dagger}Specifications\ may\ vary\ based\ upon\ configuration.$ 

 $<sup>{\</sup>it *Actual driving range will vary based on driving and battery conditions.}$ 







