

# PIPER SERIES.

The Piper series, the best selling walker in its class, was designed for those who need a lightweight yet solid and stable platform.

It's the one of the strongest aluminum walkers on the market today and is available in four sizes. Piper Super Mini with a 17" seat height weighing just 12 lbs to the Piper Tall with a 23" seat height weighing just 13 lbs.

## FEATURES.

**Aluminum Frame** Superb stability, and style

**Seat** Solid plastic base with removable fitted waterproof cushion

**Mesh Basket** Removable

**Wheels** 8" 40 mm wide, Solid grip / Longer lasting polyurethane, Suitable for indoor & outdoor use

**Colour** Vibrant Red, Sapphire Blue, Silver, Black, Forest Green



## SPECIFICATIONS.

Super Mini	Mini	Regular	Tall
			
12 lbs 16.5" Seat Width 17" Seat Height 21" Width x 26" Length 29" - 33" Handle Height 275 lbs Weight Capacity	13 lbs 16" Seat Width 19" Seat Height 21" Width x 26" Length 30" - 34" Handle Height 275 lbs Weight Capacity	13 lbs 16" Seat Width 21" Seat Height 21" Width x 26" Length 31" - 36" Handle Height 275 lbs Weight Capacity	13 lbs 16" Seat Width 23" Seat Height 23" Width x 26" Length 35" - 39" Handle Height 275 lbs Weight Capacity

## OPTIONAL ACCESSORIES.

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Mono Brakes
- Tote Bag (fits inside basket)
- Tray - Seat Mounted



EVOLUTION TECHNOLOGIES INC.  
2500 Dapple Avenue  
Port Coquitlam, BC V3C 2J9

☎ +1 800 955 2556  
info@evolution.ca  
www.evolution.ca

EVOLUTION TECHNOLOGIES USA INC.  
2997 Delta Tang Road  
Ferdale WA 95840

☎ +1 509 836 2556  
info@evolutionusa.com  
www.evolutionusa.com