

rollz

Available in four colours



Ice-Blue Deep Purple Orange Cool Grey

Accessories



Rollz Motion Cane Holder

Take your cane with you wherever you go. Because of the clever design nearly every cane fits.



Rollz Motion Cup Holder

the cup holder can be placed into position quickly, enabling you to take a water bottle along.



Rollz Motion Shopper

You can put your shopping directly in the Shopper and take them with you.



Rollz Motion Back Support

Want to take a quick rest? Use the back support for a comfortable break.



Rollz Motion Travel Cover

This accessory helps you keep the Rollz Motion including transport chair pack safe in the back of your car or even in the cargo hold of an airplane.

rollz



Triumph Mobility, Inc.

151 Bentley Ave. Unit 4
Nepean, Ontario, K2E 6T7

Tel: 1-855-546-0711
Fax: 1-866-225-3059

email: info@triumphmobility.com

triumphmobility.com

Authorized Dealer

Limited lifetime warranty.

Please visit triumphmobility.com
for additional information.



@Triumphmobility /Triumphmobility +Triumphmobility1

rollz®
takes you further



Rollator and
Transport Chair
in one



Out and about with the Rollz Motion

Going out and about without a care in the world with the Rollz Motion. Walking with the rollator or sitting in the transport chair, enjoying a lovely day in a city: Strolling through the old town, resting in the transport chair while enjoying the lush park, and then walking again to that nice restaurant to end this interesting and inspiring day.



With a flick of the wrist

With a flick of the wrist you can transform this rollator into a transport chair and back. Walk steadily, upright and comfortably, or rest along the way in the wheelchair if you feel like it.



Stable

Walking tall, steadily and supported between the ergonomic hand grips that are adjustable in height.

Comfortable

Comfortable ride and natural suspension thanks to the cushioned tires.

Easy to manoeuvre

Light driving comfort and easy manoeuvrability with the nimble front wheels.

Compact

The Rollz Motion is easily foldable into a compact package.

Ease

Resting on a comfortable seat cushion.

Assistant

Cross doorsteps and take curbs easily with the standard curb climber.

Safe

Innovative drum brakes for safe and balanced braking.

Specifications

Weight Rollator: 23.7 lb/10.6 kg

Transport Chair Pack: 8.16 lb/3.7 kg

Weight capacity: 275 lb/125 kg

Unfolded h×w×l: 34-39×25.5×24.8"
/86-102×65.5×63 cm

Width Folded: 11.5"/29 cm

Seat rollator: 16.5×8"/42 cm

Seat Transport chair: 16.5×15"/42×38 cm

Seat height: 21.5"/55 cm

Front wheel diameter: 8"/20.3 cm

Rear wheel diameter: 12"/30.5 cm

Distance between rear wheels: 22.5"/57 cm

Model Numbers

Orange: 510-2010RM0001

Ice Blue: 510-2010RM0002

Cool Grey: 510-2010RM0004

Deep Purple: 510-2010RM0006

