

XPRESSO SERIES.

The Xpresso, along with all Evolution walkers, was designed to be the strongest and most stable of all the center folding walkers on the market today.

With cable-free brakes, a center folding seat, durable folding basket, the Xpresso folds easily and compact.

FEATURES.

Innovative Design Roman arch provides support & stability

Seat Center folding durable waterproof polyurethane cushion

Brake System Cable-free / Easy to adjust

Collapsible Basket Foldable & removable

Wheels 8" 40 mm wide, Solid grip / Longer lasting polyurethane

Colors Ruby Red, Azure Blue, Chrome Silver, Ebony Black, Kivi Metallic Green, Pearlescent White, Solar Metallic Orange, Mellow Yellow, Pearlescent Champagne, Cappuccino Brown



SPECIFICATIONS.

Mini



18 lbs
18" Seat Width
16" Seat Height
22.5" Width x 26" Length
30.5" - 33.5" Handle Height
300 lbs Weight Capacity

Regular



18 lbs
18" Seat Width
21" Seat Height
22.5" Width x 26" Length
30.5" - 35.5" Handle Height
300 lbs Weight Capacity

Tall



18 lbs
18" Seat Width
24" Seat Height
22.5" Width x 26" Length
33.5" - 40" Handle Height
350 lbs Weight Capacity

OPTIONAL ACCESSORIES.

- Cup Holder
- Cane Holder
- Curb Climber
- Slow Down Brakes
- Deluxe Back Strap Pad
- Oxygen Tank Holder
- Wire Basket



EVOLUTION TECHNOLOGIES INC.
2550 Drexel Avenue
Port Coquitlam, BC V3C 2J9

☎ +1 800 955 2756
info@evolution.ca
www.evolution.ca

EVOLUTION TECHNOLOGIES USA INC.
2597 Duffin Hwy Road
Ferdale WA 95840

☎ +1 800 835 2756
info@evolutionusa.com
www.evolutionusa.com